

Bausch Arti-Fol[®] metallic Shimstock-Film - 12 microns



Arti-Fol[®]12μ metallic is a high-tech test film with distinctly improved features. This test film is made of metallic polyester film (Shimstock-film) only 12μ thick. The combination of a newly created color coating and a metallic film offers certain advantages in some applications.

Item:	Width:	Color:	Order No.:
20 m in dispenser	one-sided 22 mm	black	BA.BK30
20 m in dispenser	one-sided 22 mm	red	BA.BK31
20 m in dispenser	one-sided 22 mm	green	BA.BK32
20 m in dispenser	one-sided 22 mm	blue	BA.BK33
100 sheets (8mm x 60mm)	one-sided 8 mm	red	BA.BK35
100 sheets (8mm x 60mm)	uncoated 8 mm	-	BA.BK38
20 m in dispenser	uncoated 16 mm	-	BA.BK39
20 m	one-sided 75 mm	black	BA.BK730
20 m	one-sided 75 mm	red	BA.BK731

Bausch Arti-Fol[®] Articulating-Film ultra-thin - 8 microns



Bausch Arti-Fol[®] are especially suitable for representing static and dynamic occlusion in different colors. The first step is to inspect the concentric contact (static occlusion) in red and the second step to inspect the eccentric contacts (dynamic occlusion) in black. The color sequence can always be alternated. Several different colors can be used for a much more precise representation of dynamic occlusion.

Item:	Contents:	Width:	Color:	Order No.:	Refill-Box:
20 m in dispenser	one-sided	22 mm	black	BA.BK 20	BA.BK 1020
20 m in dispenser	one-sided	22 mm	red	BA.BK 21	BA.BK 1021
20 m in dispenser	one-sided	22 mm	green	BA.BK 22	BA.BK 1022
20 m in dispenser	one-sided	22 mm	blue	BA.BK 23	BA.BK 1023
20 m in dispenser	two-sided	22 mm	black	BA.BK 24	BA.BK 1024
20 m in dispenser	two-sided	22 mm	red	BA.BK 25	BA.BK 1025
20 m in dispenser	two-sided	22 mm	green	BA.BK 26	BA.BK 1026
20 m in dispenser	two-sided	22 mm	blue	BA.BK 27	BA.BK 1027
20 m in dispenser	one-sided	22 mm	white	BA.BK 29	BA.BK 1029

Bausch PROGRESS 100[®] 100 microns with progressive color transfer



This smooth fiber reinforced paper with high coloring capacity adapts perfectly to the occlusal surfaces. The occlusal contact points or centric contacts are marked very precisely because of the paper's progressive color transfer.

Item:	Contents:	Color:	Order No.:
Plastic dispenser	300 strips	blue	BA.BK51
Plastic dispenser	300 strips	red	BA.BK52
Plastic dispenser horseshoe	50 sheets	blue	BA.BK53
Plastic dispenser horseshoe	50 sheets	red	BA.BK54
Box with pre-cut strips	100 strips	blue	BA.BK57
Box with pre-cut strips	100 strips	red	BA.BK58

Bausch 200μ Articulating Paper with progressive color transfer



**THE LIGHTER THE BITE, THE LIGHTER THE MARK
THE HARDER THE BITE, THE DARKER THE MARK**

Item:	Contents:	Color:	Order No.:
Plastic dispenser	300 strips	blue	BA.BK01
Refill-Box	300 strips	blue	BA.BK 1001
Plastic dispenser	300 strips	red	BA.BK02
Refill-Box	300 strips	red	BA.BK 1002
Plastic dispenser horseshoe	50 sheets	blue	BA.BK03
Plastic dispenser horseshoe	50 sheets	red	BA.BK04
Box with booklets straight	300 strips	blue	BA.BK05

Bausch Arti-Check[®] 40μ micro-thin coated with liquid colours on both sides



Bausch Arti-Check[®]40μ micro-thin Articulating Papers are thin and tear-resistant papers which are coated with liquid colors on both sides. False or smear contacts can thus be avoided.

Item:	Contents:	Color:	Order No.:
Box with booklets	straight 200 strips	blue	BA.BK 09
Box with booklets	straight 200 strips	red	BA.BK 10
Box with sheets	100 x 70 mm 100 sheets	blue	BA.BK 11
Box with sheets	100 x 70 mm 100 sheets	red	BA.BK 12
Roll in dispenser	16 mm wide 15 m	blue	BA.BK 13
Roll in dispenser	16 mm wide 15 m	red	BA.BK 14
Roll in dispenser	22 mm wide 10 m	blue	BA.BK 15
Roll in dispenser	22 mm wide 10 m	red	BA.BK 16