

The Cavex Alginate system

Cavex alginates have long exceeded the criteria established by institutions such as ISO/ADA for this process. Cavex has set its own standard, which is reflected in the Cavex alginate branding.



Performance and specifications

	Cavex CA37	Cavex Orthotrace	Cavex ColorChange	Cavex Cream Alginate	Cavex Impressional
Contents	453 g	500 g	500 g	500 g	500 g
Number of imprints	approx. 35	approx. 35	approx. 38	approx. 35	approx. 38
Consistency	solid	solid	elastic	elastic	very elastic
Dimensional stability	2 days	2 days	9 days	2 days	2 days
Detail reproduction	50µm	50µm	25µm	5µm	25µm
Colour	pink	magenta	light pink	purple	blue
Colour change	no	no	yes	no	no
Setting types	FS & NS	Extra FS	FS	NS	FS & NS
Time in mouth	NS = 90 / FS = 60	60	60	90	NS = 90 / FS = 60
Total setting time	NS = 210 / FS = 150	130	150	210	NS = 210 / FS = 150

Hygiene Product



Cavex ImpreSafe starterkit:

- 1 Bottle of ImpreSafe disinfection concentrate (1 Litre),
- 1 ImpreSafe Container
- 1 Timer

Cavex ImpreSafe

A disinfecting system for alginate, silicone and polyether impressions. Kills bacteria, viruses and fungi in only 3 minutes. 100% safe for all impression surfaces. 1 litre of Cavex ImpreSafe equals 33 litres ready-to-use disinfection fluid

Container dimensions: L 30.9 x W 15.4 x H 9.9 cm



Cavex GreenClean (Alginate Remover)

Efficient alginate and gypsum cleaner; pH-neutral and 100% biodegradable. Fast, effective and harmless to metal trays. Fresh lemon fragrance. Indications: for removing alginate and gypsum residues from impression trays, spatulas and other instruments

Contents: A jar of 1 kg powder with measuring device (=500 clean trays)