

## Obtura III Max® System



### BENEFITS OF THE OBTURA III MAX®:

- Easily fills root canals with warm gutta percha in seconds for MAXimum obturation success
- Control the speed of gutta percha extrusion with the ergonomic handpiece design for superior tactile sensation
- Cartridge-Free operation resulting in low cost per canal

### PRODUCT SPECIFICATIONS:

- Thermal protector designed for patient safety
- Lightweight, ergonomic handpiece design reduces hand and wrist fatigue
- Unmatched 2 year warranty

<b>OB.823-800</b>	Obtura III Max® Heated Gutta Percha System
<b>OB.823-801</b>	Obtura Accessory Kit
<b>OB.823-813</b>	Obtura III Max® Multi-tool
<b>OB.823-815</b>	Obtura III Max® Thermal Protectors
<b>OB.823-618</b>	Obtura Plunger Seal Assembly
<b>OB.822-609</b>	Obtura Cleaning Kit

## Obtura MaxPack



### BENEFITS OF THE OBTURA MAXPACK HEAT SOURCE:

- Convenient, cordless design allows any operator to become your endodontic operator saving you time and reducing your needs to purchase multiple pieces of equipment
- Create an endodontic post space, remove excess gutta percha during lateral condensation, or compact gutta percha to prepare for warm gutta percha injection techniques, all with one device
- Confidence in quality. Backed by Obtura Spartan Endodontics comprehensive one-year warranty

### PRODUCT SPECIFICATIONS:

- Four temperature pre-sets ranging from 150°C to 230°C
- Lightweight: weighing only 2.22 ounces
- Six slots and hexagonal grip allowing heat pluggers to be inserted in any directions
- Safety: Auto shut-off after 10 minutes of downtime
- Ergonomic: cordless design allows for unlimited mobility
- High quality Lithium Ion battery minimize downtime for charging



<b>OB.824-700</b>	Obtura Spartan MaxPack System
<b>OB.824-702</b>	MaxPack Battery
<b>OB.824-719</b>	Obtura MaxPack Spreader (Narrow)- 25/.03
<b>OB.824-720</b>	Obtura MaxPack Spreader (Standard)- 25/.05
<b>OB.824-722</b>	Obtura MaxPack Plugging Stand Black (Holds 6)
<b>OB.824-723</b>	Obtura MaxPack Plugging Protector Sleeve (Bag of 10)

### OB.824-722

Obtura MaxPack Plugging Stand Black (Holds 6)