



- Incorporates a new, 35 minute 'Quick Cycle' for faster processing of instruments
- Viewing window and illuminated chamber enables the operator to make a simple check to ensure the machine is operating correctly
- Fully compliant to EN 15883 and HTM2030 - the recommended 'best practice' standard – to give you real peace of mind
- Seamless hands free process from the dentist, through cleaning to sterilization in the autoclave, eliminating sharps injuries and instrument damage

## UltraClean 3

Bench Top Washer Disinfector with Dryer

Part No.	88002110
Name	Standard Bench Top
Dimensions	
Width	560mm
Height	555mm
Depth	575mm
Max net weight	43kg
Gross weight	46kg
Chamber width	400mm
Chamber depth	450mm
Chamber height	300mm
Load capacity	
Chamber volume	60 litre
Basket volume	50 litre
Process monitoring	
Fault alert	
Warranty	
Parts (unlimited usage, excludes consumables)	1 year

Part No.	88002110
Name	Standard Bench Top
Connections	
Drain connection	25mm dia
Voltage	230 VAC
Wattage	2750 W
Frequency (hz)	50 Hz
Current	12A
Power pump	150 W
Hydraulic connection cold water	
min/max pressure	2.5 Bar min 5 Bar max
min/max temperature	+5 to +15 degrees C
Hydraulic connection warm water	
min/max pressure	2.5 Bar min 5 Bar max
min/max temperature	+5 to +15 degrees C
maximum capacity	8 kg
Noise	
Max room temperature permitted	52 dB(A) 40 degrees C

