

No.1 Rated Fiber Reinforcement Ribbon

Ribbond prevents fracture failures in dental composites and acrylics. Ribbond's unique combination of ultra-high strength fibers, enhanced bondability and patented cross-link lock-stitch leno weave makes the strength and fracture toughness of Ribbond reinforced prostheses unsurpassed by other fiber reinforcements.

Indications for Use

Periodontal Splinting (Anterior & Posterior)
 Diastema Closure Retention, Bridges (Extra-Coronal & Intra-Coronal)
 Endodontic Restorations
 Large Composite Restorations
 Denture Reinforcement And Repair
 Split Tooth Syndrome Treatment
 Over-Denture Component Reinforcement
 Orthodontic Retainers
 Traum Stbilization

Each kit comes with enough material to make about 18 to 20 canine-to-canine periodontal splints, 25-30 single pontic anterior bridges, or 40-70 fiber reinforced composite restorations. Ribbond reinforcement fibers have an indefinite shelf life.

Original-Ribbond 0.35 mm thick - Introduced in 1991, Original-Ribbond set standard in fiber reinforcement for ease of use and durability.

Ribbond-THM 0.18 mm thick - Introduced in 2001, Ribbond THM is an excellent all-purpose fiber. THM- thinner, higher modulus.

Ribbond-Ultra 0.12 mm thick - Introduced in in 2013 Ribbond ULTRA is our premium and thinnest dental fiber of all Ribbond fiber ribbons.

Ultra - ultra thin, higher modulus.



RIBBOND ASSORTED STARTER KITS - 22 CM LONG PIECES OF RIBBOND IN 2, 3 AND 4 MM SIZES

RI.AST-22 Original-Ribbond 0.35 mm thick Assorted Starter Kits
RI.ASTHM-22 Ribbond-THM 0.18 mm thick Assorted Starter Kits

- Assorted starter kits contain of 2, 3 and 4 mm sizes in 22 cm
- Special autoclavable Ribbond scissors
- Easy to understand instructions
- Dead soft tinfoil for pre-measuring in the mouth

RIBBOND ASSORTED REFILL KITS - 22 CM LONG PIECES OF RIBBOND IN 2, 3 AND 4 MM SIZES

RI.REAST-22 Original-Ribbond 0.35 mm thick Assorted Refill Kits
RI.REASTHM Ribbond-THM 0.18 mm thick Assorted Refill Kits

- Assorted starter kits contain of 2, 3 and 4 mm sizes in 22 cm
- Easy to understand instructions
- Dead soft tinfoil for pre-measuring in the mouth

RIBBOND REFILL KITS 68 CM - LONG PIECE IN TAPE

Original-Ribbond 0.35 mm thick Refill Kits

RI.RE2-68: 2 mm size
 RI.RE3-68: 3 mm size
 RI.RE4-68: 4 mm size

Ribbond-THM 0.18 mm thick Refill Kits

RI.REO7THM-68: 1 mm size (ortho)
 RI.RE2THM-68: 2 mm size
 RI.RE3THM-68: 3 mm size
 RI.RE4THM-68: 4 mm size
 RI.RE7THM-68: 7 mm size

- Each kits contain one size in 68 cm long piece
- Easy to understand instructions
- Dead soft tinfoil for pre-measuring in the mouth

RIBBOND REFILL KITS 22 CM - 3 x 22 CM LONG PIECES OF RIBBOND

Original-Ribbond 0.35 mm thick Refill Kits

RI.RE2-22: 3 x 2 mm size
 RI.RE3-22: 3 x 3 mm size
 RI.RE4-22: 3 x 4 mm size

Ribbond-THM 0.18 mm thick Refill Kits

RI.THMRE2-22: 3 x 2 mm size
 RI.THMRE3-22: 3 x 3 mm size
 RI.THMRE4-22: 3 x 4 mm size

- Each kits contain of 3 pieces in same width size of Ribbond
- Easy to understand instructions
- Dead soft tinfoil for pre-measuring in the mouth

SPECIAL AUTOCLAVABLE RIBBOND SCISSORS

RI.SCISSORS



Ribbond consists of high-strength polyethylene fibers that are comparable with those used e.g. for the manufacture of bullet-proof vests. These fibers are so sturdy that they can not be cut using a normal pair of scissors, but must instead be cut with specially honed Ribbond scissors.