

DIE STONE

ULTRA HARD

ISO TYPE: 4

DIE STONE with superior quality



Pink

ULTRA HARD **SNOW ROCK** DIE STONE UD™



DK.UD18P (Pink)

18kg Carton Box (1,5kg with Aluminum foil pack × 6pcs × 2 middle boxes)

DENTAL STONE

VERY HARD

ISO TYPE: 3

Premium DENTAL STONE with excellent properties



Yellow

HARD **SNOW ROCK** DENTAL STONE ND™



DK.ND18Y (Yellow)

18kg Carton Box (1,5kg with Aluminum foil pack × 6pcs × 2 middle boxes)

DENTAL PLASTER

STUDY & MOUNTING USE, HARD

ISO TYPE: 2

Production of study model



White

ECONOMIC **SNOW ROCK** DENTAL PLASTER™



DK.DP25

25kg Paper Bag (Plastic pack in a Paper bag)

FLOW ROCK

FLUID, Low-Expansion Base Stone

ISO TYPE: 4

Liquid pour base stone with long working time and low expansion



Yellow



Sky Blue



Pink

FLUID **SNOW ROCK** FLOW ROCK™



DK.UDF25LY (Yellow)

DK.UDF25LB (Sky Blue)

DK.UDF25LP (Pink)

11,25kg Carton Box- BULK Packing type (Plastic pack in a carton box)

ORTHODONTIC ROCK

FOR ORTHODONTIC USE & VERY HARD

ISO TYPE: 3

Esthetics and diagnosis with good workability



Pure White

ORTHODONTIC **SNOW ROCK** ORTHODONTIC ROCK™



DK.NDOT12L

11,25kg Carton Box- BULK Packing type (Plastic pack in a carton box)