

## Ultra Hydrophilic Impression Materials: Genie VPS® (Vinyl polysiloxane)

Genie® is available in a full range of wash and tray materials - five viscosity and two set time that perform like magic.



- Ultra Hydrophilic chemistry provides better results in a wet oral environment
- Two set times: 2:10 for Rapid Set, 4:30 for Standard Set
- Exceeds ISO 4380 standards for consistency, strain in compression, and detail reproduction
- Excellent flow characteristics under pressure for better surface coverage
- Exceptional batch-to-batch consistency for predictable results every time
- Clean and easy dispensing from Type 50 dispenser gun
- Berry flavored to make patients happy

Genie VPS® Impression Materials				
	Extra Light Body	Light Body	Regular Body	Heavy Body
<b>Rapid Set (2:10)</b>	SU.77600	SU.77610	SU.77620	SU.77630
<b>Standard Set (4:30)</b>	SU.77605	SU.77615	SU.77625	SU.77635

Genie VPS® includes: 2- pack Berry Flavored (2 x 50ml Cartridges, 6 Mixing Tips - XLB and LB varieties include 6 Intraoral Tips)



### Genie Putty Impression Materials

Genie® is a berry-scented, addition silicone impression material system intended for all impression techniques. The material is available in putty, heavy, regular, light, and extra-light viscosities. Both a rapid set (2:10 minutes) and a standard set (4:30 minutes) are available in each viscosity. Fabrication of gypsum casts should be delayed for 30 minutes and epoxy casts for 60 minutes to avoid porosity. Impressions may be stored for up to 2 weeks before pouring.



Rapid Set (2:30)	Standard Set (4:30)
SU.77650	SU.77655

Berry Flavored (2 x 300ml jars and 2 scoops)

### Genie® Dispenser Gun Type S-50

Genie® Dispenser Gun is used to dispense 50 ml cartridges of Genie® bite registration and VPS impression material with a ratio of 1:1.

SU.70051

## Impression Materials Genie® Magic Mix



Genie® Magic Mix is a fast, easy and economical way to use Genie® Impression Material. The new disposable 380ml cartridges have 5% more material than foil bag systems and can be used with automix dispensing machines.

	Rapid Set (2:30)	Standard Set (4:30)
<b>Regular Body</b>	SU.78800	SU.78805
<b>Heavy Body</b>	SU.78810	SU.78815

1 x 380ml cartridge, 10 mixing tips and 1 bayonet ring

## Bite Registration Material: Genie® Bite



Genie® Bite starts off with a light, fluffy consistency that stays where you put it without slumping. Once set, it becomes rock-hard, so it captures and keeps great detail without compressing or flexing during model mounting.

Rapid Set (60sec)	Standard Set (90sec)
SU.77645	SU.77640

Mint Flavor (2 x 50ml Cartridges, 6 Green Mixing Tips)

## Soothing Hand Hygiene: Moist SURE™



**A complete line of hand hygiene products.**

Designed specifically for infection prevention in dental practices, Moist SURE™ products are clinically proven to moisturize your hands, leaving them soft and velvety.



SU.95741	<b>Moist SURE™</b> Liquid Sanitizer 443 ml
----------	---

- Active ingredient: 63% isopropyl alcohol
- Kills MRSA and VRE in 5 sec



SU.95727	<b>Moist SURE™</b> Lotion Soap 1L
SU.95728	<b>Moist SURE™</b> Lotion Soap 3.8L

- Long-lasting moisturization
- Kills MRSA and VRE in 15 sec

## Personal Protective Equipment: com-fit® Face Masks



**Filtration, breathability and comfort.**

A simple, complete line of facial protection meets the needs of any office.



**com-Fit™ Super Sensitive Face Mask Ear Loop White**

- Box of 50 **SU.20346**

- Latex-free with no inks, dyes, perfumes or scratchy surfaces to irritate skin
- Provides 99% filtration at 0.1 µ



**com-Fit™ Super High Filtration Face Mask Ear Loop Blue**

- Box of 40 **SU.20326**

- 2-way DELNET® film barrier provides fluid protection with excellent breathability
- Provides 99% filtration at 0.1 µ



**com-Fit™ Groovy Mask Ear Loop**

- Box of 40 **SU.20350**

- Unique design your staff and your patients will love
- Provides 99.6% filtration at 0.1 µ



**com-Fit™ Fluid Resistant Procedural Mask Ear Loop Blue**

- Box of 50 **SU.20324**

- Contains a fluid resistant outer layer
- Provides 99% filtration at 0.1 µ

## Instrument Sterilization: Sporox® II



**High-level disinfection without irritating glutaraldehydes.**

Sporox's powerful hydrogen peroxide-based is formulated for heat-sensitive instruments that cannot be autoclaved.



No mixing,  
ready-to-use

**SU.75156 Sporox® II Sterilizing and Disinfecting Solution 3.8L**

- Ready to use – no mixing, heating or activation needed
- Safe, with no noxious odors
- Unlike glutaraldehydes, won't bond infectious proteins
- Re-use the solution up to 21 day

## Instrument Sterilization: PeelPro Self Seal Pouches



**The benefits of a high-quality sterilization pouch, but at a value price.**

Peel Pro offers solid construction and reliable performance.

SU.88000	70 x 255 mm	Box of 200
SU.88005	83 x 165 mm	
SU.88010	89 x 255 mm	
SU.88020	133 x 279 mm	
SU.88025	190 x 330 mm	

- Dual internal and external steam sterilization indicators
- Wide seal for excellent puncture protection
- Medical-grade paper allows steam to penetrate

## Instrument Cleaning Solutions: PUREVAC® Evacuation System Cleaner



**An amazing solution for both medical and dental practices.**  
The ZymeX family of enzymatic cleaners is safe and effective in processing all types of medical and dental devices.



1:125  
concentrate



No mixing,  
ready-to-use



SU.21381	<b>ZymeX™</b> Dual Enzyme Cleaner Concentrate 950 ml
SU.21380	<b>ZymeX™</b> Dual Enzyme Cleaner Concentrate 3.8L

- Applicable for clean & disinfect dental instruments
- Breaks down blood, rust and mineral deposits
- Attacks protein, carbohydrates and lipid containing soils
- Free rinsing, leaves no residue
- Ideal for ultrasonic, reprocessing and washer/sterilizers

SU.21389	<b>ZymeX™</b> Super Foaming Enzymatic Spray 660 ml
----------	--

- Ready-to-use, direct spray to container
- Long-lasting foam remains active for over 4 hours
- Dual enzyme formula breaks down blood and protein soils
- Safe for plastic, rubber, stainless steel, carbon steel and soft metals

## Evacuation System Cleaning: PUREVAC® Evacuation System Cleaner



**Powerful cleaners... made for any evacuation system!**  
Dissolves away the toughest debris to provide optimal suction. Super-concentrated Purevac® cleans your entire system—including the trap and leaves it smelling lemon fresh.



1:32  
concentrate



SU.EX21112	<b>PUREVAC®</b> Evacuation System Cleaner 2L
SU.EX21114	<b>PUREVAC®</b> Evacuation System Cleaner 5L

- Applicable for clean, disinfect & deodorant suction system
- Remove calcium and deposits that can hinder your system's performance
- Fresh lemon fragrance. Non-foaming and biodegradable

SU.21116	<b>PUREVAC®</b> Auto Walkabout Dispenser
----------	--

- Creates a turbulent mixture of air and Purevac to clean your evacuation system with ease
- Holds enough solution to clean three operatories and one cuspidor